

25mm Vapor Cell Specification

Description: 25mm Pyrex Vapor Cell


Manufacturer: Triad Technology, Inc., Longmont, CO

Fill Material :

Rubidium	(75mm standard length)
Cesium	(75mm standard length)
Potassium	(75mm standard length)
Sodium	(75mm standard length)
Iodine	(100mm standard length)

Window Material : Schott Borefloat

Dimensions :



Body Specification :

Fill Stem Height:	< 10.0mm
Thickness:	~ 1.5mm
Diameter:	25.0mm
Length:	Customer Spec.

Window Specification :

Diameter:	25.0mm
Thickness:	3.0mm
Window Angle:	0°
Wedge Angle:	N/A

Document: <http://data.sacher-laser.com/GasCell/25mmGasCell.pdf>
 Note: Specification are subject to change without further notice

